IEM SINERGIE bigh care

High tech, high care

Cutting-edge solutions for Nuclear Medicine, Radiopharmacy and

Tema Sinergie Profile

A world leading company


for radioprotection systems, pharmaceutical and health care equipment certified to quality standards since 2000.



Our History

Established in 1985, Tema Sinergie has become the world leader in designing and manufacturing of shielded isolators for Nuclear Medicine applications





Our Mission





Facility - As Is



Expanding the business

From 7.500 m² (81.000 sf) of existing offices, production area and warehouse...



Facility - As Is





Facility - To Be





Tema Sinergie Customers





Over 2000 customers Over 8000 installations

In the world.

Business Units



Nuclear Medicine

Production



Barrier Isolation Technology

Production



Radiotherapy and Medical Physics

Distribution



Nuclear Medicine





Tema Sinergie represents

the most qualified company worldwide in the production of manipulation and dispensing systems of radiopharmaceuticals

WITH MORE THAN 35 YEARS OF EXPERIENCE

Nuclear Medicine




dispensing hot cells, to dose drawing systems, radiopharmaceutical injectors, monitoring systems and much more, Tema Sinergie is the ideal partner for each and every application within Nuclear Medicine.



Customized Solutions for RP





NM References



























KAR₁₀₀ | Introduction



▼ KARℓ100 is a mobile unit designed to handle both the dispensing and injections of FDG and other PET radiopharmaceuticals using syringes and our unique automatic injection system Rad-

Inject.



KAR₁₀₀ | Introduction



KARe₁₀₀: a valuable assistant in PET







 $\mathit{KAR\ell}_{100}$ The multi-dose administration system

RAD-INJECT
The smallest automated injector for radiopharmaceuticals

Class 2b Medical Device

STAND for RAD-INJECT

KAR₁₀₀ | Workflow

3

Injection



2

Measurement, dilution and dispensing




















MULTI-TRACER COMPATIBLE

¹⁸F ⁶⁴Cu ¹²³I ⁹⁹Tc ⁶⁸Ga ¹⁷⁷Lu ⁹⁰Y ²²³Ra ¹³¹I ⁸⁹Zr ¹¹C



ASTER: the solution for RLT monodose vial injection to patient



ASTER is the safest and easiest solution to manage the administration of radiopharmaceuticals in Radiometabolic Therapy

ASTER






The system allows to automatically manage, from a single dose vial, the infusion of radiopharmaceuticals used in theranostics, minimizing exposure of the operator and avoiding any type of contamination.


Overview ASTER Set Connection





Overview Safety Transfer





SETTING

INFUSION

POST INFUSION

Pros



Only one small hole in the Glass Lead Pot for a single needle. The needle will never get in contact with the Lead Glass Pot.

No chance to get punctured by the needle.

No radiopharmaceutical leaks from the set during the disconnection from the vial.

MINIMIZE RADIO ACTIVE EXPOSURE

STERILITY

PROTECT THE OPERATOR





CRP Compact



for dispensing and sealing both open and closed vials of radiopharmaceuticals in sterile conditions. It allows safe and rapid radiopharmaceutical dispensing into vials minimizing the operator's tasks.

CRP Compact



- Automatic filling unit for open and closed vials for radiopharmaceutical production
- The most precise and reliable vial dispenser in the market



T Fully validated both for diagnostic and therapeutic radiopharmaceutical s in GMP productions


New platform release - CRP Compact +



■ Specifically design to work also with RayDyLyo® push-fit caps (or similar ones like DAIKYO PLASCAP ®)



T Designed to work also with conventional stoppers and aluminum crimp for open vial dispensing



AGLTS 2 is the band new fully automatic GMP compliant glove integrity testing system, developed referring to the Positive Pressure Decay Method which follows the international standard ISO 14644-7 Annex E.5.

A new generation device way beyond the mere execution of glove testing.

AGLTS 2



AGLTS



AGLTS



T Easy-to-use, requires no external piping for power supply, and compressed air

T Powered by the brand new easy-to replace Lithium battery

Arr Performances down to **100** $m \mu m!$



Data Transmission





T Wireless Wi-Fi data transmission based on a TCPIP protocol

The recognition of the single glove to be tested is done by means of Radio Frequency IDentification (RFID) Technology, standard feature of the equipment.



Software Compliancy



Good Automated Manufacturing Practice





Electronic records and Electronic signatures

Mr.ALESSANDRO TRERE'

alessandro.trere@temasinergie.com



TEMA CONTACT FOR THAILAND / SOUTH EAST ASIA



